



o v e s c o
innovation in scope



COAG DISSECTOR

Monopolar electrocautery instrument for precise coagulation and blunt dissection of tissue

- Concise rotation of instrument tip for optimized orientation
- Blunt tissue dissection as used in laparoscopic surgery
- Effective HF hemostasis similar to using surgical tweezers
- Small geometry for easy maneuverability





COAG DISSECTOR

The Ovesco Coag Dissector combines ease of use and safe blunt tissue dissection with precise coagulation of bleeding.


The tissue can be targeted precisely due to the easy rotation of the instrument tip.




Monopolar electro-surgical instrument for blunt tissue dissection and coagulation

- Rotatable, flexible instrument shaft for a precise alignment of the grasper tip
- Curved grasper for optimal targeting of the tissue
- Compatible with a working channel diameter of ≥ 2.8 mm
- Length of instrument shaft: 165 cm
- Variety of applications: ESD, Hemostasis, POEM,...


ESD – blunt tissue dissection



Marking of lesion




Injection to lift up lesion




Blunt tissue dissection by opening Coag Dissector


Hemostasis



Grasping of bleeding vessel



HF cauterization



Effective hemostasis

Ovesco Endoscopy AG
Dorfackerstr. 26
72074 Tuebingen/Germany

Phone +49(0)7071 96528-160
Fax +49(0)7071 96528-260
service@ovesco.com

Ovesco Endoscopy USA Inc.
120 Quade Dr
Cary, NC 27513
USA

Phone +1 919 651 9449
Fax +1 408 608 2077
customerservice@ovesco-usa.com



Ovesco Endoscopy AG is a medical device company specializing in the fields of flexible endoscopy and endoluminal surgery. Ovesco develops, manufactures, and markets innovative products for the treatment of gastrointestinal disease. Ovesco products stand for therapeutic efficacy and efficient application.

www.ovesco.com