


HIGH LOW  
**PROTECTION LEVEL**



**BLUE**

**GREEN**

**PINK**

**YELLOW**

**WHITE**

**WELDING POINT**

·Ultrasonically welded to the mask, makes the soft knit earloops more comfortable and secure

**FLAT EARLOOP**

·Flat earloop provides comfort even wear for a long time

**INNER LAYER**


·Soft, special designed material to reduce skin irritation and allergic reactions  
 ·All days comfort provided

**FILTER**

·Highest quality filter media for best ASTM rating,filtration value and breathability

**OUTER LAYER**

·Added protection with unique material for maximum fluid resistance / protection against body fluids



**TOPSAFE EARLOOP FACE MASK**

Meets all current ASTM 2100 performance criteria.

CLASSIFICATION	TYPE IIR EN14683	BFE(%)	PFE(%)	DELTA P (H <sub>2</sub> O/CM <sup>2</sup> )	MM Hg <sup>☆</sup>	FLAME SPREAD
<b>TOPSAFE ELITE</b>	○	≥99.8	≥99.3	< 5.0	160	Class I
<b>TOPSAFE PLUS</b>	○	≥98.0	≥99.3	< 5.0	160	Class I
<b>DISMED PREMIER</b>	△	≥95.0	≥98.9	< 4.0	80	Class I

○ TYPE IIR EN14683 certified    △ None    ☆ Resistance to penetration by synthetic blood

Earloop Face Mask Ref. Codes		
COLOR		TOPSAFE ELITE
BLUE		1A3100
GREEN		1A3101
YELLOW		1A3102
PINK		1A3103
WHITE		1A3104



Packing: 50pcs/box, 40boxes/ctn  
 Dimension/ctn: 52cm\*40cm\*38.5cm  
 Gross Weight/ctn: 9.2kg

