

DISPOSABLE ENDOSCOPIC DISTAL ATTACHMENT

Distal Cap

A visualisation cap, used on the distal end of an endoscope and available in a range of sizes to fit most diagnostic and therapeutic endoscopes. The soft material protects the tissue while introducing the endoscope and maintains a good field of view. They are widely used in Colonoscopy procedures and help to unfold the colon, plus EMR, ESD and POEM surgeries.

Advantages:

- Maintains a good field of view during endoscopy.
- Helps to unfold the colon during colonoscopy procedures.
- Provides more working space in surgery.
- Protects the tissue.
- Compatible with a range of scope sizes.

